

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08613 - Goat Anti-ABCA2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ABCA2, ATP-binding cassette, sub-family A (ABC1), member 2, ABC2, MGC129761, ATP-binding cassette, sub-family A, member 2, OTTHUMP00000064733

Official Symbol: ABCA2

Accession Number(s): NP_001597.2; NP_997698.1

Human GeneID(s): [20](#)

Non-Human GeneID(s): 11305 (mouse), 79248 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_001597.2; NP_997698.1).

Immunogen

Peptide with sequence C-KKQSDNLEQQETEP, from the internal region of the protein sequence according to NP_001597.2; NP_997698.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:2000.

Western blot: Preliminary experiments gave an approx 35kDa band in Human Brain (Amygdala, Frontal Cortex, Hippocampus and Substantia Nigra) lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 270kDa according to NP_001597.2 and of 273kDa according to .NP_997698.1. The 35kDa band was successfully blocked by incubation with the immunizing peptide.

IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 3.75µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat

EB08613 (3.75µg/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.