

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05674-B - Goat Anti-ACADM, Biotinylated Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ACADM, acyl-CoA dehydrogenase, C-4 to C-12 straight chain, ACAD1, MCAD, MCADH, acyl-Coenzyme A dehydrogenase, C-4 to C-12 straight chain, FLJ18227, FLJ93013, FLJ99884, medium-chain acyl-CoA dehydrogenase, medium-chain specific acyl-CoA dehydrogenase

Official Symbol: ACADM

Accession Number(s): NP_000007.1; NP_001120800.1; NP_001272971.1;
NP_001272972.1; NP_001272973.1

Human GeneID(s): [34](#)

Important Comments: This antibody is expected to recognise both reported isoforms (NP_000007.1; NP_001120800.1).

Immunogen

Peptide with sequence C-RLIVAREHIDKYKN., from the C Terminus of the protein sequence according to NP_000007.1; NP_001120800.1; NP_001272971.1; NP_001272972.1; NP_001272973.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 45kDa band observed in Human Heart lysates (calculated MW of 46.6kDa according to NP_000007.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 2-4µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



Biotinylated EB05674 (3µg/ml) staining of Human Heart lysate (35µg protein in RIPA buffer), exactly mirroring its parental non-biotinylated product. Primary incubation was 1 hour. Detected by chemiluminescence, using streptavidin-HRP and using NAP blocker as a substitute for skimmed milk.