

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB10204 - Goat Anti-ACHRE / CHRNE Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** AchR epsilon subunit, ACHRE, cholinergic receptor, nicotinic, epsilon, cholinergic receptor, nicotinic, epsilon polypeptide, CMS1D, CMS1E, CMS2A, FCCMS, nicotinic acetylcholine receptor epsilon polypeptide, SCCMS, CHRNE

**Official Symbol:** CHRNE

**Accession Number(s):** NP\_000071.1

**Human GeneID(s):** [1145](#)

**Non-Human GeneID(s):** 11448 (mouse), 29422 (rat)

### Immunogen

Peptide with sequence C-DQEATGEEVSD, from the Internal region of the protein sequence according to NP\_000071.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested


**Peptide ELISA:** antibody detection limit dilution 1:1000.

**Western blot:** Approx 55kDa band observed in Mouse Fetal Skeletal Muscle lysates (calculated MW of 54.9kDa according to NP\_033733.1). Recommended concentration: 1-3µg/ml.

### Species Reactivity

**Tested:** Mouse

**Expected from sequence similarity:** Human, Mouse, Rat, Cow



EB10204 (1µg/ml) staining of Mouse Fetal Skeletal Muscle lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.