

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08017 - Goat Anti-ADRB3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ADRB3, adrenergic, beta-3-, receptor, BETA3AR, Beta-3 Adrenergic Receptor

Official Symbol: ADRB3

Accession Number(s): NP_000016.1

Human GeneID(s): [155](#)

Immunogen

Peptide with sequence PHENSSLAPWPD-C, from the N Terminus of the protein sequence according to NP_000016.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 45kDa band observed in Human Brain (Cortex) lysates (calculated MW of 43.5kDa according to NP_000016.1). Recommended concentration: 0.03-0.1µg/ml. Primary incubation was 1 hour.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB08017 (0.03µg/ml) staining of Human Brain (Cerebral Cortex) lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.