

## UK Office

**Everest Biotech Ltd**  
Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07943 - Goat Anti-BHMT Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** BHMT, betaine-homocysteine methyltransferase

**Official Symbol:** BHMT

**Accession Number(s):** NP\_001704.1

**Human GeneID(s):** [635](#)

**Non-Human GeneID(s):** 12116 (mouse), 81508 (rat)

### Immunogen

Peptide with sequence C-EQQLKELFEKQK, from the C Terminus of the protein sequence according to NP\_001704.1 .

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:16000.

**Western blot:** Approx 45kDa band observed in Rat Liver lysates (calculated MW of 44.9kDa according to human NP\_001704.1, 45.0kDa according to mouse NP\_057877.1 and rat NP\_110477.1). This product has been successfully used in WB on Mouse (PMID:28605831, 23807810 and 22209966). In transfected HEK293 transiently expressing BHMT a band of approx. 49kDa is observed. This band is not observed in the non-transfected HEK293. This product has been successfully used in WB on Human (PMID: 27320863). Recommended concentration: 0.03-0.1µg/ml. Primary incubation was 1 hour.

**IHC:** In paraffin embedded Human Kidney shows strong cytoplasm staining in PCT cells. Recommended concentration: 2-4µg/ml.

### Species Reactivity

**Tested:** Human, Mouse, Rat

**Expected from sequence similarity:** Human, Mouse, Rat, Dog

### Specific References

**This antibody has been successfully used in Western blot on Mouse:**

Jeon J, Oh JJ, Kwak HC, Yun H, Kim HC, Kim YM, Oh SJ, Kim SK.  
Age-Related Changes in Sulfur Amino Acid Metabolism in Male C57BL/6 Mice.  
Biomol Ther (Seoul). 2018 Mar 1;26(2):167-174  
PMID: 28605831

**This antibody has been successfully used in Western blot on Human and Rat:**

Ahn CW, Jun DS, Na JD, Choi YJ, Kim YC.  
Alleviation of hepatic fat accumulation by betaine involves reduction of homocysteine via up-regulation of methylenetetrahydrofolate dehydrogenase.  
Biochem Biophys Res Commun. 2016 Aug 26;477(3):440-7.  
PMID: 27320863

**This antibody has been successfully used on Rat:**

Kim YC, Kwon do Y, Kim JH.

Alterations in the metabolomics of sulfur-containing substances in rat kidney by betaine.  
Amino Acids. 2014 Apr;46(4):963-8.

PMID: 24390397

**This antibody has been successfully used in Western blot on Mouse:**

Kwon DY, Jung YS, Kim SJ, Kim YS, Choi DW, Kim YC.

Alterations in sulfur amino acid metabolism in mice treated with silymarin:  
a novel mechanism of its action involved in enhancement of  
the antioxidant defense in liver.

Planta Med. 2013 Aug;79(12):997-1002.

PMID: 23807810

**This antibody has been successfully used on Mouse:**

Yun KU, Ryu CS, Lee JY, Noh JR, Lee CH, Lee HS, Kang JS, Park SK, Kim BH, Kim SK.

Hepatic metabolism of sulfur amino acids in db/db mice.

Food Chem Toxicol. 2013 Mar;53:180-6.

PMID: 23220616

**This antibody has been successfully used in Western blot on Mouse:**

Yun KU, Ryu CS, Oh JM, Kim CH, Lee KS, Lee CH, Lee HS, Kim BH, Kim SK.

Plasma homocysteine level and hepatic sulfur amino acid metabolism in mice fed a  
high-fat diet.

Eur J Nutr. 2012 Jan 1.

PMID: 22209966

**This antibody has been successfully used on Human:**

Ryu CS, Kwak HC, Lee KS, Kang KW, Oh SJ, Lee KH, Kim HM, Ma JY, Kim SK.

Sulfur amino acid metabolism in doxorubicin-resistant breast cancer cells.

Toxicol Appl Pharmacol. 2011 Jun 14.

PMID: 21703291

EB07943 (0.03µg/ml) staining of Rat Liver lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.  
HEK293 overexpressing BHMT RC203148) and probed with EB07943 (mock transfection in first lane), tested by  
Origene.