

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB12680 - Goat Anti-CENPE (aa2363-2374) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CENPE, centromere protein E, 312kDa, CENP-E, KIF10, PPP1R61, Centromere autoantigen E (312kD), centromere-associated protein E, kinesin family member 10, kinesin-related protein CENPE, protein phosphatase 1, regulatory subunit 61

**Official Symbol:** CENPE

**Accession Number(s):** NP\_001804.2; NP\_001273663.1

**Human GeneID(s):** [1062](#)

**Important Comments:** This antibody is expected to recognize both reported isoforms (NP\_001804.2; NP\_001273663.1).

### Immunogen

Peptide with sequence C-TQDNKNPHVTSR, from the internal region of the protein sequence according to NP\_001804.2; NP\_001273663.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:128000.

**Western blot:** Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 316kDa according to NP\_001804.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human