

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06105 - Goat Anti-ARAP3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CENTD3, centaurin, delta 3, ARAP3, FLJ21065, ARF-GAP, RHO-GAP, ankyrin repeat and plekstrin homology domains-containing protein 3, PtdIns(3,4,5)P3-binding protein, phosphoinositide binding protein, ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3, DRAG1, Arf and Rho GAP adapter protein 3

Official Symbol: CENTD3

Accession Number(s): NP_071926.4

Human GeneID(s): [64411](#)

Immunogen

Peptide with sequence CTSSPPSSQPLT, from the C Terminus of the protein sequence according to NP_071926.4.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: This antibody has been successfully used in WB on scirrhous human gastric cancer cell lines.

Immunoprecipitation: This antibody has been successfully used in IP on scirrhous human gastric cancer cell lines.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Dog, Pig, Cow