

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08494 - Goat Anti-CHRNA1 / ACHRB Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CHRNA1, ACHRB, cholinergic receptor, nicotinic, beta 1 (muscle), CHRNA, CMS1D, CMS2A, SCCMS, cholinergic receptor, nicotinic, beta polypeptide 1 (muscle), nicotinic acetylcholine receptor beta 1 subunit

Official Symbol: CHRNA1

Accession Number(s): NP_000738.2

Human GeneID(s): [1140](#)

Immunogen

Peptide with sequence C-QEQEDHDALKED, from the internal region of the protein sequence according to NP_000738.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.


Western blot: Approx 55kDa band observed in human brain (cerebellum, hippocampus, frontal cortex) lysates (calculated MW of 56.7kDa according to NP_000738.2).

Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Cow



EB08494 (1 μ g/ml) staining of human cerebellum lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.