



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08389 - Goat Anti-Dachshund homolog 2 / DACH2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DACH2, dachshund homolog 2 (Drosophila), RP11-345E19.1, FLJ31391, MGC138545, dachshund 2

Official Symbol: DACH2

Accession Number(s): NP_444511.1

Human GeneID(s): [117154](#)

Non-Human GeneID(s): 93837 (mouse)

Immunogen

Peptide with sequence C-TRKQAVNSSRPGR, from the internal region of the protein sequence according to NP_444511.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Approx 65-70kDa band observed in lysates of cell line NIH3T3 (calculated MW of 65.3kDa according to human NP_444511.1 and of 68.6kDa according to mouse NP_291083.1). Recommended concentration: 0.03-0.1µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse

EB08389 (0.03µg/ml) staining of NIH3T3 cell lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.