

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05072 - Goat Anti-DAPP1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DAPP1, dual adaptor of phosphotyrosine and 3-phosphoinositides, BAM32, DKFZp667E0716, HSPC066, B-cell adapter molecule of 32 kDa, b lymphocyte adapter protein Bam32, dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide, hDAPP1

Official Symbol: DAPP1

Accession Number(s): AAF14578.1

Human GeneID(s): [27071](#)

Important Comments: This antibody was raised against the C-terminus of NP_055210.1. However, NCBI has changed the version to NP_055210.2 thus losing the C-terminal epitope.

Immunogen

Peptide with sequence SALCISPEEKTDHK, from the C Terminus of the protein sequence according to AAF14578.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Approx 30kDa band observed in Human Kidney lysates (calculated MW of 32.0kDa according to AAF14578.1). Recommended concentration: 0.5-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB05072 (0.5µg/ml) staining of Human Kidney lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.