

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10886 - Goat Anti-DBF4B (aa193-207) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DBF4B, DBF4 homolog B (S. cerevisiae), ASKL1, DRF1, FLJ13087, MGC15009, ZDBF1B, chifb, ASK-like protein 1, Dbf4-related factor 1, activator of S phase kinase-like protein 1, activator of S-phase kinase-like protein 1, chifon homolog B, protein DBF4 homolog B, zinc finger, DBF-type containing 1B

Official Symbol: DBF4B

Accession Number(s): NP_663696.1; NP_079380.1

Human GeneID(s): [80174](#)

Important Comments: This antibody is expected to recognize both reported isoforms.

Immunogen

Peptide with sequence QLSLASLCVKKQPK, from the internal region of the protein sequence according to NP_663696.1; NP_079380.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human