

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB09705 - Goat Anti-DTX2 / Deltex-2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** DTX2, deltex homolog 2 (Drosophila), KIAA1528, MGC71098, RNF58, OTTHUMP00000024490, deltex 2, ring finger protein 58, zinc ion binding protein

**Official Symbol:** DTX2

**Accession Number(s):** NP\_065943.2; NP\_001096066.1

**Human GeneID(s):** [113878](#)

**Non-Human GeneID(s):** 304591 (rat)

**Important Comments:** This antibody is expected to recognize both reported isoforms (NP\_065943.2; NP\_001096066.1). Reported variants NP\_001096065.1, NP\_065943.2, NP\_001096064.1 represent identical protein.

### Immunogen

Peptide with sequence C-EKLSTASGYSDVTD, from the internal region of the protein sequence according to NP\_065943.2; NP\_001096066.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.


**Western blot:** Preliminary experiments gave an approx 85kDa band in lysates of cell lines A549, HepG2 and NIH3T3 after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 67.2kDa according to NP\_065943.2. The 85kDa band was successfully blocked by incubation with the immunizing peptide.

**IHC:** Paraffin embedded Human Small Intestine. Recommended concentration: 5µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Rat, Cow



EB09705 (5µg/ml) staining of paraffin embedded Human Small Intestine. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.