

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05321 - Goat Anti-DUSP10 / MKP5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DUSP10, MKP5, dual specificity phosphatase 10, MKP-5, dual specificity phosphatase MKP-5, serine/threonine specific protein phosphatase

Official Symbol: DUSP10

Accession Number(s): NP_009138.1; NP_653329.1; NP_653330.1

Human GeneID(s): [11221](#)

Important Comments: This antibody is expected to recognise both reported Human isoforms. The variants NP_653329.1 and NP_653330.1 represent identical protein.

Immunogen

Peptide with sequence CRILTPKLMGVETVV, from the C Terminus of the protein sequence according to NP_009138.1; NP_653329.1; NP_653330.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 53kDa band seen in Human Liver lysate. Recommended for use at 0.5-2µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB05321 staining (0.5µg/ml) of Human Liver lysate (RIPA buffer, 30µg total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.