

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05636 - Goat Anti-EAR2 / NR2F6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: NR2F6, nuclear receptor subfamily 2, group F, member 6, EAR2, EAR-2, ERBAL2, ERBA-related gene-2, v-erb-a avian erythroblastic leukemia viral oncogene homolog-like 2

Official Symbol: NR2F6

Accession Number(s): NP_005225.2

Human GeneID(s): [2063](#)

Non-Human GeneID(s): 13864 (mouse), 245980 (rat)

Immunogen

Peptide with sequence C-STFNWPYGSGQ, from the C Terminus of the protein sequence according to NP_005225.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Consistent band observed in Human Testis at approx 55-60kDa [Predicted MW of approx. 45kDa according to NP_005225] with 1 hour incubation of primary.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow