

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06592 - Goat Anti-ERK2 / MAPK1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MAPK1, ERK2, ERK, p38, p40, p41, ERT1, MAPK2, PRKM1, PRKM2, P42MAPK, p41mapk, mitogen-activated protein kinase 1, protein tyrosine kinase ERK2, mitogen-activated protein kinase 2, extracellular signal-regulated kinase 2

Official Symbol: MAPK1

Accession Number(s): NP_002736.3; NP_620407.1

Human GeneID(s): [5594](#)

Non-Human GeneID(s): 26413 (mouse)

Important Comments: This antibody is expected to recognise an epitope corresponding to aa 9-19 of both reported protein isoforms of human ERK2 / MAPK1 protein.

Immunogen

Peptide with sequence CAAGPEMVRGQVF, from the internal region of the protein sequence according to NP_002736.3; NP_620407.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 38kDa band observed in HepG2 lysates (calculated MW of 41.4kDa according to NP_002736 and NP_620407). Recommended concentration: 0.5-2µg/ml. Primary incubation was 1 hour.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow

EB06592 (0.5µg/ml) staining of HepG2 lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.