

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06116 - Goat Anti-FAZF / TZFP / ZBTB32 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ZBTB32, TZFP, FAZF, testis zinc finger protein, FAXF, zinc finger and BTB domain containing 32, Rog, ZNF538, OTTHUMP00000045790, repressor of GATA

Official Symbol: ZBTB32

Accession Number(s): NP_055198.1

Human GeneID(s): [27033](#)

Immunogen

Peptide with sequence RPSTSPCCPSSSTT, from the C Terminus of the protein sequence according to NP_055198.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: A faint approx 55kDa band observed in Human Bone Marrow lysates (calculated MW of 53.0kDa according to NP_055198). Recommended concentration: 2-5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB06116 (2 μ g/ml) staining of Human Bone Marrow lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.