

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06421 - Goat Anti-FOXP4 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** FOXP4, forkhead box P4, winged-helix repressor FOXP4, fork head-related protein like A, FLJ40908, FLJ44184, hFKHLA, OTTHUMP00000043536

**Official Symbol:** FOXP4

**Accession Number(s):** NP\_001012426.1; NP\_612466.1; NP\_001012427.1

**Human GeneID(s):** [116113](#)

**Non-Human GeneID(s):** 74123 (mouse)

**Important Comments:** This antibody is expected to recognize all three reported isoforms (NP\_001012426.1; NP\_612466.1; NP\_001012427.1).

### Immunogen

Peptide with sequence VESASETIRSAPSG-C, from the N Terminus of the protein sequence according to NP\_001012426.1; NP\_612466.1; NP\_001012427.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:4000.

**Western blot:** Customer's experiments gave a band at approx 75kDa in lysates of FOXP4-transfected HEK293 (calculated MW of 73.5kDa according to NP\_001012426.1). Recommended concentration: 1-3µg/ml.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Cow