

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10777 - Goat Anti-galectin-1 (mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Lgals1, lectin, galactose binding, soluble 1, AA410090, Gal-1, Galbp, L-14.5, L14, Lect14, galectin-1, S-Lac lectin 1, beta-galactoside binding protein, beta-galactoside-binding lectin L-14-I, galaptin, galectin-1, lactose-binding lectin 1

Official Symbol: Lgals1

Accession Number(s): NP_032521.1

Human GeneID(s): [3956](#)

Non-Human GeneID(s): 16852 (mouse), 56646 (rat)

Immunogen

Peptide with sequence KPGECLKVRGEVASD, from the N Terminus (near) of the protein sequence according to NP_032521.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.


Western blot: Approx 15kDa band observed in lysates of Mouse and Rat Pancreas and of cell line NIH3T3 (calculated MW of 14.9kDa according to NP_032521.1).

Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Pig, Cow



EB10777 (0.1 μ g/ml) staining of Mouse (A) and Rat (B) Pancreas lysates and NIH3T3 (C) lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.