

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07564 - Goat Anti-GALR2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GALR2, galanin receptor 2, GALNR2, MGC125983, MGC125984

Official Symbol: GALR2

Accession Number(s): NP_003848.1

Human GeneID(s): [8811](#)

Non-Human GeneID(s): 14428 (mouse), 29234 (rat)

Immunogen

Peptide with sequence C-RYPLHSRELTPRN, from the internal region of the protein sequence according to NP_003848.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Approx 40kDa band observed in Human Heart lysates (calculated MW of 41.7kDa according to NP_003848.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB07564 (0.3µg/ml) staining of Human Heart lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.