

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08246 - Goat Anti-GLIPR1 / RTVP-1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GLIPR1, GLI pathogenesis-related 1 (glioma), RTVP-1, CRISP7, GLIPR, RTVP1, glioma pathogenesis-related protein, glioma pathogenesis-related protein 1, related to testis-specific, vespid, and pathogenesis proteins 1, testes-specific vespid and pathogenesis protei

Official Symbol: GLIPR1

Accession Number(s): NP_006842.2

Human GeneID(s): [11010](#)

Immunogen

Peptide with sequence C-NRQRDQVKRYYS, from the internal region of the protein sequence according to NP_006842.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary experiments gave an approx 40kDa band in human prostate lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 30.4kDa according to NP_006842.2. The ~40kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human