

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07945 - Goat Anti-GPR3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GPR3, G protein-coupled receptor 3, ACCA, OTTHUMP00000043209, adenylate cyclase constitutive activator

Official Symbol: GPR3

Accession Number(s): NP_005272.1

Human GeneID(s): [2827](#)

Non-Human GeneID(s): 14748 (mouse), 266769 (rat)

Important Comments: This antibody is not expected to cross-react with GPR6 and GPR12.

Immunogen

Peptide with sequence C-LQRHLLPASH, from the internal region of the protein sequence according to NP_005272.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 35kDa band observed in Human Brain (Hippocampus and Frontal Cortex) lysates (calculated MW of 35.0kDa according to NP_005272.1). Recommended concentration: 0.3-1µg/ml. An additional band of unknown identity was also consistently observed at 30kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog



EB07945 (0.3µg/ml) staining of Human Brain (Frontal Cortex) lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.