

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11808 - Goat Anti-histamine H3 receptor Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: G protein-coupled receptor 97, GPCR97, G-protein coupled receptor 97, H3R, HH3R, histamine H3 receptor, histamine receptor H3, HRH3, OTTHUMP00000031459

Official Symbol: HRH3

Accession Number(s): NP_009163.2

Human GeneID(s): [11255](#)

Non-Human GeneID(s): 99296 (mouse), 85268 (rat)

Immunogen

Peptide with sequence C-SVTRAVSYRAQQGDT, from the internal region of the protein sequence according to NP_009163.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 48kDa band observed in lysates of cell line A549 (calculated MW of 48.7kDa according to NP_009163.2). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB11808 (0.1 μ g/ml) staining of A549 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.