

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12202 - Goat Anti-HOXB6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HOXB6, homeobox B6, HOX2, HOX2B, HU-2, Hox-2.2, homeo box 2B, homeo box B6, homeobox protein Hox-2.2, homeobox protein Hox-2B, homeobox protein Hox-B6, homeobox protein Hu-2

Official Symbol: HOXB6

Accession Number(s): NP_061825.2

Human GeneID(s): [3216](#)

Non-Human GeneID(s): 15414 (mouse), 497986 (rat)

Immunogen

Peptide with sequence EPRKSDCAQDKS, from the internal region of the protein sequence according to NP_061825.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.


Western blot: Approx 26kDa band observed in lysates of cell lines HEK293, CACO-2 and MRC5 PD19 (calculated MW of 25.4kDa according to NP_061825.2). Recommended concentration: 1-3µg/ml.

IHC: In paraffin embedded Human Lung shows nuclear staining in the alvioli. Recommended concentration, 5-10µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB12202 (5 μ g/ml) staining of paraffin embedded Human Lung. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



EB12202 (1 μ g/ml) staining of HEK293 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.