

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11909 - Goat Anti-HOXD13 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HOXD13, homeobox D13, BDE, BDSD, HOX4I, SPD, homeo box 4I, homeo box D13, homeobox protein Hox-4I, homeobox protein Hox-D13

Official Symbol: HOXD13

Accession Number(s): NP_000514.2

Human GeneID(s): [3239](#)

Non-Human GeneID(s): 15433 (mouse), 288154 (rat)

Immunogen

Peptide with sequence CTKDQPQGSHFWKS, from the internal region of the protein sequence according to NP_000514.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Approx 38kDa band observed in Human (Cerebellum), Mouse and Rat Brain lysates (calculated MW of 36.1kDa according to Human NP_000514.2, 35.9kDa according to Mouse NP_032301.2 and 34.9kDa according to Rat NP_001099356.1). Recommended concentration: 0.1-1µg/ml. Primary incubation 1 hour at room temperature.

Species Reactivity

Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Pig, Cow

EB11909 (0.3µg/ml) staining of Human (A), Mouse (B) and (0.1ug/ml) Rat (C) Brain lysate (35µg protein in RIPA buffer). Detected by cemiluminescence.