

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07673 - Goat Anti-KLK2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: KLK2, kallikrein-related peptidase 2, KLK2A2, MGC12201, hK2, glandular kallikrein 2, kallikrein 2, prostatic

Official Symbol: KLK2

Accession Number(s): NP_005542.1; NP_001002231.1

Human GeneID(s): [3817](#)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_005542.1; NP_001002231.1).

Immunogen

Peptide with sequence C-KHQSLRPDEDSSHD, from the internal region of the protein sequence according to NP_005542.1; NP_001002231.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 22kDa band observed in Human Prostate lysates (calculated MW of 24.7kDa according to NP_001002231.1). Recommended concentration: 0.1-0.3µg/ml.

IHC: In paraffin embedded Human Prostate shows vesicular staining in the secretory cells. Recommended concentration: 2.5-3.8µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB07673 (0.1µg/ml) staining of Human Prostate lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.

EB07673 (3.8µg/ml) staining of paraffin embedded Human Prostate. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.