

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05965 - Goat Anti-MDA5 / IFIH1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** IFIH1, MDA5, MDA-5, Hlcd, MGC133047, IDDM19, melanoma differentiation associated protein-5, DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide, interferon induced with helicase C domain 1, helicard

**Official Symbol:** IFIH1

**Accession Number(s):** NP\_071451.2

**Human GeneID(s):** [64135](#)

### Immunogen

Peptide with sequence SNGYSTDENFRYL-C, from the N Terminus of the protein sequence according to NP\_071451.2.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**IHC:** In paraffin embedded Human Pancreas, shows strong cytoplasmic staining of some cells in the islets of Langerhans. Recommended concentration: 10ug/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human

### Specific Reference

**This antibody has been successfully used in WB in Human:**

Sato S, Hoshino K, Satoh T, Fujita T, Kawakami Y, Fujita T, Kuwana M.

RNA helicase encoded by melanoma differentiation-associated gene 5 is a major autoantigen in patients with clinically amyopathic dermatomyositis: Association with rapidly progressive interstitial lung disease.

Arthritis Rheum. 2009 Jul;60(7):2193-200.

PMID: 19565506

EB05965 (10µg/ml) staining of paraffin embedded Human Pancreas. Microwaved antigen retrieval with Tris/EDTA buffer pH9, HRP-staining.