

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05619 - Goat Anti-MPS1 / RPS27 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MPS1, RPS27, ribosomal protein S27 (metallopanstimulin 1), MPS-1, metallopanstimulin 1, 40S ribosomal protein S27, RP11-422P24.3, metallopan-stimulin 1

Official Symbol: RPS27

Accession Number(s): NP_001021.1

Human GeneID(s): [6232](#)

Non-Human GeneID(s): 57294 (mouse), 94266 (rat)

Immunogen

Peptide with sequence PLAKDLLHPSPEE-C, from the N Terminus of the protein sequence according to NP_001021.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx. 12kDa band observed in lysates of Human T-cell line Jurkat (calculated MW of 9.5kDa according to NP_001021.1). Recommended concentration: 0.03-0.1µg/ml. Primary incubation was 1 hour.

IHC: Paraffin embedded Human Tonsil. Recommended concentration: 3.75µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat

EB05619 (0.03µg/ml) staining of Jurkat lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.
EB05619 (3.75µg/ml) staining of paraffin embedded Human Tonsil. Steamed antigen retrieval with citrate buffer
pH 6, AP-staining.