

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB11953 - Goat Anti-myeloperoxidase Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** MPO, myeloperoxidase

**Official Symbol:** MPO

**Accession Number(s):** NP\_000241.1

**Human GeneID(s):** [4353](#)

**Important Comments:** The immunizing peptide represents the N terminus of the mature protein.

### Immunogen

Peptide with sequence C-DKAYKERRESIKQ, from the internal region (near N terminus) of the protein sequence according to NP\_000241.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:128000.

**Western blot:** Approx 85kDa band observed in Human Peripheral Blood

Mononucleocytes lysates (calculated MW of 83.9kDa according to NP\_000241.1).

Recommended concentration: 0.3-1µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human

EB11953 (0.3µg/ml) staining of Human Peripheral Blood Lymphocyte lysate (35µg protein in RIPA buffer).  
Primary incubation was 1 hour. Detected by chemiluminescence.