

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09175 - Goat Anti-NDEL1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: NDEL1, nudE nuclear distribution gene E homolog (A. nidulans)-like 1, DKFZp451M0318, EOPA, MITAP1, NUDEL, NUDE-like protein, endooligopeptidase A, mitosis-associated protein MITAP1, nudE nuclear distribution gene E homolog like 1

Official Symbol: NDEL1

Accession Number(s): NP_001020750.1; NP_110435.1

Human GeneID(s): [81565](#)

Non-Human GeneID(s): 83431 (mouse), 170845 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_001020750.1; NP_110435.1).

Immunogen

Peptide with sequence CRNFAKDQASRKSY, from the internal region of the protein sequence according to NP_001020750.1; NP_110435.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 38+39kDa bands observed in Human Brain (Cerebellum) lysates (calculated MW of 37.0kDa according to NP_001020750.1 and of 38.4kDa according to NP_110435.1). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow

EB09175 (0.1µg/ml) staining of Human Cerebellum lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.