

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10189 - Goat Anti-NDUFS2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: complex 1, mitochondrial respiratory chain, 49-KD subunit, NADH dehydrogenase (ubiquinone) Fe-S protein 2, NADH dehydrogenase (ubiquinone) Fe-S protein 2, 49kDa (NADH-coenzyme Q reductase), NADH-ubiquinone oxidoreductase NDUFS2 subunit, NDUFS2

Official Symbol: NDUFS2

Accession Number(s): NP_004541.1; NP_001159631.1

Human GeneID(s): [4720](#)

Non-Human GeneID(s): 226646 (mouse), 289218 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_004541.1; NP_001159631.1).

Immunogen

Peptide with sequence C-SRGVRQWQPDVE, from the internal region of the protein sequence according to NP_004541.1; NP_001159631.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 48kDa band observed in Human, Mouse and Rat Heart lysates (calculated MW of 51.9kDa according to NP_001159631.1). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow, Pig

EB10189 (0.1µg/ml) staining of Human Heart lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.