

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08251 - Goat Anti-Netrin G1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: NTNG1, netrin G1, KIAA0976, Lmnt1, axon guidance molecule, laminin 1, netrin G1 isoform G1b, netrin G1 isoform G1m, netrin G1 isoform G1n, netrin G1 isoform G1o, netrin G1f

Official Symbol: NTNG1

Accession Number(s): NP_055732.2

Human GeneID(s): [22854](#)

Immunogen

Peptide with sequence C-QSATWKEYPKPLQV, from the internal region of the protein sequence according to NP_055732.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments in Human Brain (Cerebellum, Hippocampus and Cerebral Cortex) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow, Horse