

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07284 - Goat Anti-Peroxiredoxin 2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** peroxiredoxin 2, HGNC:9353, PRDX2, MGC4104, NKEFB, PRP, PRXII, TDPX1, TSA, natural killer-enhancing factor B, thiol-specific antioxidant 1, thioredoxin peroxidase 1, thioredoxin-dependent peroxide reductase 1, torin

**Official Symbol:** PRDX2

**Accession Number(s):** NP\_005800.3; NP\_859427.1

**Human GeneID(s):** [7001](#)

**Important Comments:** This antibody is expected to recognise isoforms a and b only (NP\_005800.3 and NP\_859427.1 respectively).

### Immunogen

Peptide with sequence C-NVDDSKKEYFSKHN, from the C Terminus of the protein sequence according to NP\_005800.3; NP\_859427.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:16000.

**Western blot:** Approx 20-22kDa band observed in Human Brain lysates (calculated MW of 21.9kDa according to NP\_005800.3). Recommended concentration: 0.3-1µg/ml.

Primary incubation was 1 hour.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat



EB07284 (1 $\mu$ g/ml) staining of Human Brain lysate (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.