

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08670 - Goat Anti-PKD1 / Polycystin 1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** PKD1, polycystic kidney disease 1 (autosomal dominant), PBP, polycystic kidney disease-associated protein, polycystin 1

**Official Symbol:** PKD1

**Accession Number(s):** NP\_001009944.1; NP\_000287.2

**Human GeneID(s):** [5310](#)

**Important Comments:** This antibody is expected to recognize both reported isoforms (NP\_001009944.1 and NP\_000287.2).

### Immunogen

Peptide with sequence C-RTPLRAKNKVHP, from the C Terminus of the protein sequence according to NP\_001009944.1; NP\_000287.2.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:64000.

**IHC:** In paraffin embedded Human Kidney shows textured cytoplasmic staining towards the basal membrane in PCT. Recommended concentration: 3-5µg/ml. Paraffin embedded Human Spleen. Recommended concentration: 3.75µg/ml.

**Immunoprecipitation:** This antibody was successfully used on Human (PMID:25574838, and PMID: 29706351). **Immunofluorescence:** This antibody was successfully used in IF on Human (Durand et al, Sci Rep. 2016 Oct 24;6:35963, PMID: 27775029).

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human

### Specific References

**This antibody has been successfully used in IP on Human:**

Cornec-Le Gall E, Olson RJ, Besse W, Heyer CM, Gainullin VG, Smith JM, Audrézet MP, Hopp K, Porath B, Shi B, Baheti S, Senum SR, Arroyo J, Madsen CD, Férec C, Joly D, Jouret F, Fikri-Benbrahim O, Charasse C, Coulibaly JM, Yu AS, Khalili K, Pei Y, Somlo S, Le Meur Y, Torres VE; Genkyst Study Group; HALT Progression of Polycystic Kidney Disease Group; Consortium for Radiologic Imaging Studies of Polycystic KidneyDisease, Harris PC.

Monoallelic Mutations to DNAJB11 Cause Atypical Autosomal-Dominant Polycystic Kidney Disease.

Am J Hum Genet. 2018 May 3;102(5):832-844

PMID: 29706351

**This antibody has been successfully used in IF on Human:**

Durand N, Bastea LI, Long J, Döppler H, Ling K, Storz P.

Protein Kinase D1 regulates focal adhesion dynamics and cell adhesion through Phosphatidylinositol  
Sci Rep. 2016 Oct 24;6:35963.  
PMID: 27775029

**This antibody has been successfully used in IP on Human:**

Gainullin VG, Hopp K, Ward CJ, Hommerding CJ, Harris PC.

Polycystin-1 maturation requires polycystin-2 in a dose-dependent manner.

J Clin Invest. 2015 Jan 9. pii: 76972.

PMID: 25574838

EB08670 (3.75µg/ml) staining of paraffin embedded Human Kidney. Steamed antigen retrieval with citrate buffer  
Ph 6, AP-staining.

EB08670 (3.75µg/ml) staining of paraffin embedded Human Spleen. Steamed antigen retrieval with citrate buffer  
pH 6, AP-staining.