

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07277 - Goat Anti-Psoriasin / S100A7 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: S100A7, S100 calcium binding protein A7 (psoriasin 1), HGNC:10497, PSOR1, S100A7c, S100 calcium-binding protein A7, S100 calcium-binding protein A7 (psoriasin 1), psoriasin 1, Psoriasin

Official Symbol: S100A7

Accession Number(s): NP_002954.1

Human GeneID(s): [6278](#)

Immunogen

Peptide with sequence C-HKYTRRDDKIDKP, from the internal region of the protein sequence according to NP_002954.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 15kDa band observed in Human Skin lysates (calculated MW of 11.5kDa according to NP_002954.1). Recommended concentration: 1-3µg/ml. In transfected HEK293 transiently expressing full-length Human S100A7 (myc and DYKDDDDK tagged), a band of approx. 16kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-myc tag antibody. Recommended concentration, 0.1-0.5µg/ml. Approx 15kDa band was also observed in Human Skin lysates (calculated MW of 11.5kDa according to NP_002954.1). Recommended concentration: 0.3-1µg/ml.

IHC: Paraffin embedded Human Tonsil. Recommended concentration: 3.75µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB07277 (0.1µg/ml) staining of Human Skin lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.

HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human S100A7 with C-terminal MYC tag probed with EB07277 (0.5ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with EB07277 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.

EB07277 (3.75µg/ml) staining of paraffin embedded Human Tonsil. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.