

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08311 - Goat Anti-RAB8A Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RAB8A, RAB8A, member RAS oncogene family, MEL, RAB8, mel transforming oncogene, mel transforming oncogene (RAB8 homolog), mel transforming oncogene (derived from cell line NK14), mel transforming oncogene (derived from cell line NK14)- RAB8 homolog, ras-associated protein RAB8

Official Symbol: RAB8A

Accession Number(s): NP_005361.2

Human GeneID(s): [4218](#)

Immunogen

Peptide with sequence C-NQGVKITPDQQKR, from the C Terminus of the protein sequence according to NP_005361.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 25kDa band observed in Human Brain (Cerebellum), Kidney and Heart lysates (calculated MW of 23.7kDa according to NP_005361.2). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Dog, Cow

EB08311 (0.1µg/ml) staining of Human Brain (Cerebellum) lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.