

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05686 - Goat Anti-RAD51C Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RAD51C, RAD51 homolog C (*S. cerevisiae*), RAD51L2, yeast RAD51 homolog 3, DNA repair protein RAD51 homolog 3, MGC104277, RAD51 homolog C, RAD51 homolog C, isoform 1

Official Symbol: RAD51C

Accession Number(s): NP_002867.1; NP_478123.1

Human GeneID(s): [5889](#)

Important Comments: This antibody is expected to recognise all reported human isoforms.

Immunogen

Peptide with sequence RGKTFRFEMQRDL-C, from the N Terminus of the protein sequence according to NP_002867.1; NP_478123.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 45kDa band seen in HeLa lysate [Predicted MW of 41kDa according to NP_478123]. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB05686 staining (2µg/ml) of HeLa lysate (RIPA buffer, 30µg total protein per lane). Primary incubated for 12 hour. Detected by western blot using chemiluminescence.