

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05587 - Goat Anti-RNF31 / ZIBRA (C terminus) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RNF31, ZIBRA, FLJ10111, FLJ23501, zinc in-between-ring-finger ubiquitin-associated domain, ring finger protein 31, MGC19975, zinc in-between-RING-finger ubiquitin-associated domain protein

Official Symbol: RNF31

Accession Number(s): NP_060469.4

Human GeneID(s): [55072](#)

Immunogen

Peptide with sequence C-EEVPLGQSIPRRRK, from the C Terminus of the protein sequence according to NP_060469.4.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 125kDa band observed in lysates of cell line JurkatT(calculated MW of 120kDa according to AAH12077.1). Recommended concentration: 0.5-1.5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Dog, Pig

EB05587 (0.5µg/ml) staining of JURKAT lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.