

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08369 - Goat Anti-RPE65 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: mRPE65, sRPE65, retinal pigment epithelium-specific protein 65kDa, LCA2, RP20, rd12, RBP-binding membrane protein, p63, retinal pigment epithelium-specific 61 kDa protein, retinal pigment epithelium-specific protein (65kD), retinitis pigmentosa 20

Official Symbol: RPE65

Accession Number(s): NP_000320.1

Human GeneID(s): [6121](#)

Immunogen

Peptide with sequence C-EVKKNARKAPQPE, from the internal region of the protein sequence according to NP_000320.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 65kDa band observed in Rat Retina lysates (calculated MW of 60.9kDa according to Human NP_000320.1 and 61.0kDa according to Rat NP_446014.1). Recommended concentration: 0.5-1.5µg/ml. Primary incubation was 1 hour.

Species Reactivity

Tested: Rat

Expected from sequence similarity: Human, Rat, Dog

EB08369 (0.1µg/ml) staining of Rat Retina lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.