

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08028 - Goat Anti-RPL17 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RPL17, ribosomal protein L17, MGC117162, rpL23, 60S ribosomal protein L17, gene encoding putative NFkB activating protein

Official Symbol: RPL17

Accession Number(s): NP_000976.1; NP_001030178.1

Human GeneID(s): [6139](#)

Non-Human GeneID(s): 319195 (mouse), 291434 (rat)

Important Comments: Both reported variants represent identical protein (NP_000976.1; NP_001030178.1)

Immunogen

Peptide with sequence RYSLDPENPTKSC, from the N Terminus of the protein sequence according to NP_000976.1; NP_001030178.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 23kDa band observed in Human Prostate lysates (calculated MW of 21.4kDa according to NP_000976.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog

Specific Reference

Smolock EM, Korshunov VA, Glazko G, Qiu X, Gerloff J, Berk BC.

Ribosomal protein L17, RpL17, is an inhibitor of vascular smooth muscle growth and carotid intima formation.

Circulation. 2012 Nov 13;126(20):2418-27.

PMID: 23065385



EB08028 (0.3µg/ml) staining of Human Prostate lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.