

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07474 - Goat Anti-SLC22A16 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** SLC22A16, solute carrier family 22 (organic cation transporter), member 16, RP1-261K5.1, CT2, FLIPT2, OCT6, OKB1, dJ261K5.1, OTTHUMP00000040478, WUGSC:RG331P03.1, carnitine transporter 2, fly-like putative organic ion transporter 2, organic cation transporter 6, solute carrier family 22, member 16

**Official Symbol:** SLC22A16

**Accession Number(s):** NP\_149116.2

**Human GeneID(s):** [85413](#)

### Immunogen

Peptide with sequence CSRNKRENTSSL, from the internal region of the protein sequence according to NP\_149116.2.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested


**Peptide ELISA:** antibody detection limit dilution 1:128000.

**IHC:** In paraffin embedded Human testis shows membranous staining in the seminiferous tubules. Recommended concentration, 1-2µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human



EB07474 (4μg/ml) staining of paraffin embedded Human Testis. Microwaved antigen retrieval with Tris/EDTA buffer pH9, HRP-staining.