

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07026 - Goat Anti-SMO (C Terminus) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SMO, smoothened homolog (Drosophila), HGNC:11119, Gx, SMOH, seven transmembrane helix receptor, smoothened

Official Symbol: SMO

Accession Number(s): NP_005622.1

Human GeneID(s): [6608](#)

Immunogen

Peptide with sequence C-NLMDTELMADSD, from the C Terminus of the protein sequence according to NP_005622.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 60kDa band observed in Human Bone Marrow lysates (calculated MW of 86.4kDa according to NP_005622.1). The observed molecular weight corresponds to findings with EB07444, a product of different design reacting to a different part of the same protein. Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Dog, Cow



A) EB07026 and B) EB07444 (0.3µg/ml) staining of Human Bone Marrow lysate (35µg protein in RIPA buffer).
Primary incubation was 1 hour. Detected by chemiluminescence.