

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07213 - Goat Anti-STAM2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** STAM2, signal transducing adaptor molecule (SH3 domain and ITAM motif) 2, HGNC:11358, DKFZp564C047, Hbp, STAM2A, STAM2B, STAM-like protein containing SH3 and ITAM domains 2, signal transducing adaptor molecule 2

**Official Symbol:** STAM2

**Accession Number(s):** NP\_005834.3

**Human GeneID(s):** [10254](#)

**Non-Human GeneID(s):** 56324 (mouse)

**Important Comments:** No cross-reaction expected with STAM1.

### Immunogen

Peptide with sequence C-QSIDPTDSKPDSQD, from the internal region of the protein sequence according to NP\_005834.3.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:8000.

**Western blot:** Approx 70kDa band observed in Human Spleen lysates (calculated MW of 58.2kDa according to NP\_005834.3). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Pandey et al J Biol Chem. 2000 Dec 8;275(49):38633-9; PMID: 10993906). Recommended concentration: 1-3µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Dog



EB07213 (1 $\mu$ g/ml) staining of Human Spleen lysate (35 $\mu$ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.