

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09861 - Goat Anti-STAT4 (aa 613-627) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: OTTHUMP00000163553, signal transducer and activator of transcription 4, SLEB11, STAT4

Official Symbol: STAT4

Accession Number(s): NP_003142.1

Human GeneID(s): [6775](#)

Immunogen

Peptide with sequence C-DHSESGEVRFHVSVEP, from the internal region of the protein sequence according to NP_003142.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: In transfected HEK293 transiently expressing full-length Human STAT4 (MYC and DYKDDDDK tagged), a band of approx.85kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-myc tag antibody. Recommended concentration, 0.2-0.5µg/ml.

Purchase the +ve control! The transfected lysate used in QC is available for sale. Please email sales@everestbiotech.com to request datasheet, pricing and delivery information.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Dog, Pig



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human STAT4 (RC206892) with C-terminal MYC tag probed with EB09861 (0.5ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with EB09861 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.