

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05699 - Goat Anti-TAB182 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TNKS1BP1, tankyrase 1 binding protein 1, 182kDa, TAB182, KIAA1741, tankyrase 1-binding protein of 182 kDa, FLJ45975, tankyrase 1-binding protein 1

Official Symbol: TNKS1BP1

Accession Number(s): NP_203754.2

Human GeneID(s): [85456](#)

Immunogen

Peptide with sequence KVSTLRESSAMAS-C, from the N Terminus of the protein sequence according to NP_203754.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.


Western blot: Preliminary experiments in Human Ovary and Lung lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).

IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB05699 (5µg/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.