

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09559 - Goat Anti-TFB1M Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TFB1M, transcription factor B1, mitochondrial, CGI-75, CGI75, mtTFB, mtTFB1, OTTHUMP00000017473, S-adenosylmethionine-6-N\1, N\1-adenosyl(rRNA) dimethyltransferase 1, homolog of yeast mitochondrial transcription factor B, mitochondrial 12S rRNA dimethylase 1

Official Symbol: TFB1M

Accession Number(s): NP_057104.2

Human GeneID(s): [51106](#)

Immunogen

Peptide with sequence C-ELKRRKSKNEEKE, from the C Terminus of the protein sequence according to NP_057104.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Approx 38kDa band observed in lysates of cell line A431 (calculated MW of 39.5kDa according to NP_057104.2). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Dog



EB09559 (0.3 μ g/ml) staining of A431 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.