

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05656 - Goat Anti-TFB2M Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TFB2M, transcription factor B2, mitochondrial, Hkp1, FLJ22661, FLJ23182, hypothetical protein FLJ23182

Official Symbol: TFB2M

Accession Number(s): NP_071761.1

Human GeneID(s): [64216](#)

Immunogen

Peptide with sequence CAYKWLYDETLEDR, from the C Terminus of the protein sequence according to NP_071761.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Preliminary experiments gave an approx 40kDa band in Human Brain lysate at 0.3 ug/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 46kDa according to NP_071761.

IHC: Paraffin embedded Human Colon. Recommended concentration: 5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB05656 (5µg/ml) staining of paraffin embedded Human Colon. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.