

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06323 - Goat Anti-TIP30 / CC3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CC3, TIP30, HTATIP2, HIV-1 Tat interactive protein 2, 30kDa, Tat-interacting protein (30kD), HIV-1 Tat interactive protein 2, 30 kD, HIV-1 Tat interactive protein 2, 30 kDa, FLJ26963, SDR44U1, short chain dehydrogenase/reductase family 44U, member 1

Official Symbol: HTATIP2

Accession Number(s): NP_001091990.1; NP_006401.3; NP_001091991.1; NP_001091992.1; NP_001091993.1

Human GeneID(s): [10553](#)

Important Comments: This antibody is expected to recognise all three reported isoforms. Variants (NP_006401.3; NP_001091991.1; NP_001091992.1) encode the same isoform.

Immunogen

Peptide with sequence AETEALSKLREDFR-C, from the N Terminus of the protein sequence according to NP_001091990.1; NP_006401.3; NP_001091991.1; NP_001091992.1; NP_001091993.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 28kDa band observed in lysates of cell line HeLa (calculated MW of 27.0kDa according to NP_006401.3). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human

EB06323 (0.6µg/ml) staining of HeLa lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.