

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06323 - Goat Anti-TIP30 / CC3 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CC3, TIP30, HTATIP2, HIV-1 Tat interactive protein 2, 30kDa, Tat-interacting protein (30kD), HIV-1 Tat interactive protein 2, 30 kD, HIV-1 Tat interactive protein 2, 30 kDa, FLJ26963, SDR44U1, short chain dehydrogenase/reductase family 44U, member 1

**Official Symbol:** HTATIP2

**Accession Number(s):** NP\_001091990.1; NP\_006401.3; NP\_001091991.1; NP\_001091992.1; NP\_001091993.1

**Human GeneID(s):** [10553](#)

**Important Comments:** This antibody is expected to recognise all three reported isoforms. Variants (NP\_006401.3; NP\_001091991.1; NP\_001091992.1) encode the same isoform.

### Immunogen

Peptide with sequence AETEALSKLREDFR-C, from the N Terminus of the protein sequence according to NP\_001091990.1; NP\_006401.3; NP\_001091991.1; NP\_001091992.1; NP\_001091993.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Approx 28kDa band observed in lysates of cell line HeLa (calculated MW of 27.0kDa according to NP\_006401.3). Recommended concentration: 0.3-1µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human

EB06323 (0.6µg/ml) staining of HeLa lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour.  
Detected by chemiluminescence.