

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08026 - Goat Anti-TPD52L1 / D53 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TPD52L1, hD53, tumor protein D52-like 1, D53, MGC8556, TPD52L2, tumor protein D52-like 2, tumor protein D53

Official Symbol: TPD52L1

Accession Number(s): NP_003278.1; NP_001003395.1 ; NP_001003396.1; NP_001003397.1

Human GeneID(s): [7164](#)

Important Comments: This antibody is expected to recognise isoforms 1, 3 and 4 (NP_003278.1, NP_001003396.1 and NP_001003397.1 resp.).

Immunogen

Peptide with sequence C-EPLQGTDEDAVASAD, from the internal region of the protein sequence according to NP_003278.1; NP_001003395.1 ; NP_001003396.1; NP_001003397.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 26kDa band observed in Human Breast cancer lysates (calculated MW of 22.5kDa according to NP_003278.1). Recommended concentration: 0.1-0.3µg/ml. See also our D52 and D54 products.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08026 (0.1µg/ml) staining of Human Breast cancer lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.