

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11011 - Goat Anti-TPH2 (aa16-29) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TPH2, tryptophan hydroxylase 2, ADHD7, FLJ37295, MGC138871, MGC138872, NTPH, neuronal tryptophan hydroxylase, tryptophan 5-hydroxylase 2, ADHD7, FLJ37295, MGC138871, MGC138872, neuronal tryptophan hydroxylase, NTPH, tryptophan 5-hydroxylase 2, tryptophan 5-monooxygenase 2, tryptophan hydroxylase 2, TPH2

Official Symbol: TPH2

Accession Number(s): NP_775489.2

Human GeneID(s): [121278](#)

Non-Human GeneID(s): 216343 (mouse), 317675 (rat)

Immunogen

Peptide with sequence RGFSLDSAVPEEHQ-C, from the N Terminus of the protein sequence according to NP_775489.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 60kDa band observed in lysates of cell line HEK293 (calculated MW of 56.1kDa according to NP_775489.2). Recommended concentration: 0.5-2µg/ml. Primary incubation was 1 hour.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow

EB11011 (0.5µg/ml) staining of HEK293 lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.